



**Thermo Scientific  
Class II, Type A2  
Biological Safety Cabinets**

Thermo Scientific Herasafe KS  
Thermo Scientific 1300 Series A2

containment. comfort.  
**convenience.**

# containment. comfort.

Our Thermo Scientific biological safety cabinets (BSCs) are designed to maximize sample protection and operator safety – every day.

For more than four decades, we have been a pioneer and innovator of BSCs. Our global experience has helped to redefine BSC design, performance and safety – changing the way you work – for the better.

Our BSCs are backed by our technology leadership and market vision to provide safe and reliable cabinets for a wide range of laboratory applications.

**5** YEAR  
**SmartPromise**  
**Comprehensive Warranty**  
delivers total peace of mind!

## technology leadership

Our technical advancements include:

1965

Our first line of clean-air products

1978

Our first NSF-certified BSC

1985

World's first TÜV-certified BSC

1993

World's first BSC with motorized front window and aerosol-tight sealing

# convenience.™

## global market vision

Our design philosophy is grounded in compliance with current worldwide safety standards.

## quality built, premium features

Our BSCs feature German engineering for outstanding reliability and quality. Discover why laboratories around the globe rely on Thermo Scientific BSCs every day.

Thermo Scientific Herasafe KS  
Biological Safety Cabinet

**ultimate protection, comfort and performance for advanced applications**



Thermo Scientific 1300 Series A2  
Biological Safety Cabinet

**exceptional efficiency, safety and value for daily applications**



1995

Our first  
LNE-certified BSC

1996

Our first DIN  
12980-certified  
BSC for cytotoxic  
drugs

2002

World's first BSC  
with innovative  
brushless  
dual-DC motor  
technology  
(Herasafe KS)

2007

Our launch of the  
energy efficient  
1300 Series A2  
BSC

2011

Our first BSC  
with SmartPort  
tubing organizer



Independently tested and certified to  
NSF/ANSI 49 for Class II biological  
safety cabinets

# innovations in containment

## first-to-market with brushless DC motors

In 2002, we pioneered the first innovative brushless DC motor technology. Now available in our full line of BSCs, our DC motor advancements have redefined and transformed BSC energy efficiency, safety, performance and reliability.

### Use less energy, reduce costs

DC motors use less energy to operate, so there is less heat emitted back into the lab – delivering better comfort and lower air conditioning costs. On average, our BSCs use 75% less energy than traditional AC motor cabinets.

### Exceptional reliability

Our dual-motor BSC design offers an expected lifetime of more than 100,000 hours – at least double the published life of traditional AC motor and single DC motor systems – resulting in significant savings on motor replacements over the life of the cabinet.

### Long filter life

Larger mini-pleat H14 HEPA filters and balanced airflow, coupled with eight-inch work aperture, results in higher filtration efficiency and longer filter life than traditional BSCs.




## check the watts!

**Thermo Scientific biological safety cabinets lead the industry in operational cost efficiency with energy consumption as low as 150W.**

## optimizing total cost of operation

Our technology innovations in brushless DC motors, Thermo Scientific SmartFlow airflow compensation and balanced filters have led to safer BSCs that are easier to use and significantly less costly to operate.



\* Comparison based on 15-year life of 1.2 m wide biological safety cabinets and 2011 average costs of energy and cooling in the U.S. "Traditional BSC" estimates based on published product specifications for BSCs using AC motor technology.

**The more you use it, the more you save!**


# and efficiency

advanced design for peace of mind

A BSC's primary function is to provide cleanliness and containment. Available on our full range of BSCs, our safety innovations – **SmartFlow™**, **Digital Airflow Verification** and **Night Set-Back Mode** – maximize safety and minimize sample contamination.

- ▶ H14 HEPA filters increase filtration efficiency at 99.995%
- ▶ Metal plenums for long life and lasting containment
- ▶ Negative-pressure, dual side walls protect from leaks to the exterior environment

**SmartFlow:** automated balance of RPM, voltage, and airflow



## SmartFlow ensures safety

**Unique airflow system takes safety and containment to new levels**

- Smart DC motors monitor and control fan speed in real time to maintain constant airflow during filter loading or temporary airflow.
- Exclusive to our BSCs, independent supply and exhaust blowers automatically balance both the downflow and inflow velocities to maintain both user and product protection.

## Digital Airflow Verification (DAVe) validates protection

**Monitoring system alerts user when containment is compromised**

- Independent sensors detect pressure changes across the exhaust and downflow plenums.
- Alarm signals when a 20% change in inflow/exhaust or downflow occurs.
- Airflow velocities are displayed on the control panel for monitoring and recording.

## Night Set-Back Mode saves energy, maintains sterility

**When the front sash is closed, our intelligent speed control automatically reduces blower speed.**

- Continuous airflow maintains a sterile working environment, even when the cabinet is not in use.
- Less air flows through the HEPA filters, extending filter life and decreasing energy consumption.
- Night Set-Back is significantly more energy efficient than similar features on competing biological safety cabinets.



# designed for comfort and

Thermo Scientific BSCs offer an efficient workspace, improved ergonomics, and added conveniences that provide the user greater productivity and ultimately, better laboratory results.



Sloped front window allows fatigue-free working posture  
Bright open chamber improves visibility  
Quiet airflow design offers noticeably reduced noise level

Front air-grill arm rest extends user reach inside cabinet

**New!** Thermo Scientific Ergolign  
Saddle Stool: Exclusive seat design  
reduces back strain and fatigue

- Clean-room-ready, with 0.5  $\mu\text{m}$  filtration system within the seat, trapping 99.9% of particles exhausted, fulfilling Class 100 requirements
- Fully adjustable seat height and tilt, using a patented icon plate mechanism
- Stores easily under the BSC to save floor space when not in use

# convenience

work smarter, more efficiently

- ▶ Easy to use
- ▶ Easy to clean
- ▶ Extremely quiet to operate

## smarter design



### NEW SmartPort!


Clean and safe routing for vacuum tubing and cables through the side of the BSC for improved organization and work efficiency. (1300 Series A2 only)

### SmartClean™ window design

Our unique window design easily lowers for thorough cleaning of inner surfaces and protects operator by maintaining inflow, even when the window is lowered. (1300 Series A2 only)

### SmartClean Plus

Fully open, front-hinged window design allows convenient posture for thorough, easy cleaning and disinfection, with full access to BSC interior. (Herasafe KS only)



### Worry-free decontamination

Optional UV light is programmable from 30 minutes to 24 hours, in 30-minute increments, extending bulb life and saving energy.

Powerful cross-beam UV irradiation in both side walls illuminates entire work area, ensuring thorough disinfection of the complete chamber (Herasafe KS only).

# Thermo Scientific Herasafe KS

## ultimate protection, comfort and performance

for advanced applications

Featuring an advanced design for easy operation and maximum safety, our Herasafe KS BSCs may be used wherever protection is paramount, including research, pharmaceutical quality control labs, pharmacies and clinics or hospitals.

### Containment

- **SmartFlow design** automates airflow compensation and balancing
- **Digital Airflow Verification (DAVe)** alerts user to compromised containment
- **Large H14 HEPA filters** remove airborne contaminants
- **Performance Factor safety indicator** displays performance status
- **Cross-beam UV disinfection system** illuminates entire work area without shadows
- **SampleGuard aerosol-tight window sealing** for added sample protection

### Comfort

- **Reduced noise level** enhances user focus
- **10° sloped front** for exceptional comfort and ergonomics
- **Spacious work area** with comfortable armrest
- **Transparent side-glass windows** maximize light and visibility

### Convenience

- **SmartClean Plus** hinged-front window design for easy cleaning
- **Intuitive eye-level control panel** with handheld controls for ultimate ADA compliance
- **Motorized frameless window sash** ensures safe sash height at the push of a button
- **Scratch-free, 316 stainless steel work surfaces** are available as single or segmented modules and are easy to clean or autoclave
- **Programmable cross-beam UV light** extends bulb life and saves energy



**SampleGuard**  
aerosol-tight window sealing  
provides extra containment  
when window is closed

### Efficiency

- **DC motors** offer low energy consumption and greater than 100,000 hours of use.
- **Night Set-Back Mode** reduces airflow to 30% during low-usage periods
- **Long filter life**
- **Long UV bulb life**





Cross-beam UV

Intuitive eye-level control panel

Transparent side-glass windows

Motorized sash

Optional universal piping/rear wall service taps

316 stainless steel interior

Customizable work surfaces

Optional electric adjustable height stand

## More choices.

Choice of 3-, 4-, 5- or 6-foot widths  
Available with 10-inch work aperture  
Choice of electrical or manual adjustable height stands

# Thermo Scientific 1300 Series A2

**optimized efficiency, safety and value**  
for your daily applications

Our Thermo Scientific 1300 Series A2 Class II, Type A2 biological safety cabinets combine safety, energy efficiency and value for today's everyday laboratory and clinical research applications.

## Containment

- **SmartFlow design** automates airflow compensation and balancing
- **Digital Airflow Verification (DAVe)** alerts user to compromised containment
- **Large H14 HEPA filters** remove airborne contaminants
- **SmartFlow indicator** displays performance status

## Comfort

- **Low noise level** enhances user focus
- **10° sloped front** for exceptional comfort and ergonomics
- **Spacious work area** with comfortable armrest
- **Increased light intensity** for a brighter work environment

## Convenience

- **SmartClean** front window design
- **Intuitive control panel** for at-a-glance monitoring
- **SmartPort** organizes tubing and cables
- **Programmable UV light** extends bulb life and saves energy

**Intuitive Control Panel:** Displays valuable safety and performance data with constant read-out of downflow/inflow velocities



### Airflow Read-Out & Hour Counter:

Clearly displays realtime performance data, including inflow and downflow velocities and hours of operation

### SmartFlow Indicator:

Visually demonstrates power and capacity to maintain protection

### Front Window Status:

Visual and audible alarm indicates whether front window is in correct working position

### Night Set-Back Indicator:

Displays reduced speed operation when front window is closed

### DAVe Airflow Indicator:

Audible and visual alarms indicate when airflow is safe or restricted



**Thermo Scientific Ergolign Saddle Stool**

Relieves back pressure and optimizes working posture during long work sessions

## More choices.

- Available in 3-, 4-, 5- or 6-foot widths
- Choice of 8- or 10-inch work aperture
- Choice of SmartCoat interior for glare-free workspace or seamless one-piece stainless steel
- Easy-to-order package configurations with accessories (manual adjustable height stand, UV light and armrest)
- Electrical adjustable height stands are also available

### 1300 Series A2

shown with optional manual floor stand, floor-stand hanging shelf, vacuum valve and rear-wall tubing clips

## Efficiency

- **DC motors** offer low energy consumption and more than 100,000 hours of use.
- **Night Set-Back Mode** reduces airflow to 30% during low-usage periods
- **Long filter life**
- **Long UV bulb life**


# Thermo Scientific Herasafe KS

## specifications

Size/Width			3-foot	4-foot	5-foot	6-foot
<b>Dimensions</b>	Exterior dimensions WxHxD	inches	39.4 x 62.4 x 31.5	51.2 x 62.4 x 31.5	63.0 x 62.4 x 31.5	74.8 x 62.4 x 31.5
	Interior dimensions WxHxD	inches	35.4 x 30.7 x 24.8	47.2 x 30.7 x 24.8	59.1 x 30.7 x 24.8	70.9 x 30.7 x 24.8
	Shipping dimensions <sup>1</sup> WxHxD	inches	43.7 x 69.3 x 36.4	55.5 x 69.3 x 36.4	67.3 x 69.3 x 36.4	79.1 x 69.3 x 36.4
	Working height of front window	inches	10			
	Maximum opening height of front window	inches	30.4			
<b>Weight</b>	Net weight	lbs	375	441	507	617
	Shipping weight <sup>1</sup>	lbs	419	496	573	694
<b>Electrical Data</b>	Voltage	V	120			
	Frequency	Hz	60			
<b>Filter Specification</b>			H14 HEPA EN 1822, 99.995% MPPS			
<b>Certification</b>			NSF/ANSI 49, ETL, CE			
<b>Ergonomics and Utilities</b>	Lighting power	fc	>70	>70	>110	>120
	Receptacles		2 GFI duplex, one on each sidewall			
	Service valves		Up to 4 may be installed through access ports (2 on each side) Up to 6 more may be factory-installed on rear wall (optional)			
<b>10-Inch Work Aperture</b>	Sound pressure level	dB (A)	<65			
	Energy consumption, operating set point <sup>2</sup>	W	210	275	350	435
	Energy consumption, Night Set-Back mode	W	40	40	70	70
	Heat output, operating mode (non-vented)	BTU/hr	717	938	1194	1484
	Heat output, Night Set-Back mode (lights off, non-vented)	BTU/hr	136	136	239	239
	Exhaust / Inflow air volume, recirculated	cfm	255	342	428	514
	Exhaust volume, thimble ducted	cfm	332	444	556	668

<sup>1</sup> Export packaging dimensions and weight upon request

<sup>2</sup> Clean filters, fans operating at nominal set point, interior lighting activated



# ordering info

The Herasafe KS standard configuration includes cabinet, UV light, armrest, 316 stainless steel indented work surface, and battery back-up for motorized sash. Stands must be ordered separately.


Size/Width		3-foot		4-foot		5-foot		6-foot	
		Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.
<b>Standard Configuration</b>	Includes cabinet, cross-beam UV, one-piece indented worktray and battery back-up for front window. Stands are sold separately.	51027116	13-998-000	51027117	13-998-001	51027118	13-998-002	51027120	13-998-003



Cross-beam UV light



Stainless steel armrests



One-piece indented worktray

Size/Width		3-foot		4-foot		5-foot		6-foot	
		Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.
<b>Stands: (Must be ordered separately)</b>	Manual adjustable floor stand, 30 to 38 inches in 2-inch increments	50116443	13-998-006	50109312	13-998-007	50116444	13-998-008	50109313	13-998-009
	Electric adjustable floor stand, with cable remote control, 30 to 38 inches working height, 120V	50119096	13-998-004	50119126	15-453-300	50119127	13-998-005	50119128	15-453-301
	Castor set (may not be used with electric stands; fixed 30-inch work surface height)	50109977	13-998-032	50109977	13-998-032	50109977	13-998-032	50109977	13-998-032



Electric adjustable stand



Manual stand



Castors

\*Not available in Canada.

# Thermo Scientific Herasafe KS

## customer-installed accessories

	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Description
<b>Exhaust Transitions</b>	50078796	13-998-022	Thimble duct exhaust transition for Herasafe KS 9 (3 ft) cabinet
	50073770	13-998-023	Thimble duct exhaust transition for Herasafe KS 12 (4 ft) or KS 15 (5 ft) cabinet
	50078795	13-998-024	Thimble duct exhaust transition for Herasafe KS 18 (6 ft) cabinet
	191780	191780	Alnor exhaust alarm for thimble duct, 120v (independent from BSC, plugs into wall outlet)
<b>User Convenience Equipment</b>	50078900	13-998-027	Utility basket for use with hanging bars 50073943 or 1911312
	50078901	13-998-028	2 Utility baskets including additional hanging bar
	50073943	13-998-029	Hanging bar with set of 6 hooks for KS 9 (3 ft) cabinet only
	1911312	15-474-212	Hanging bar with set of 12 hooks for KS 12, 15, 18 (4, 5, 6 ft.) cabinets
	1911313	15-474-213	Replacement stainless steel armrest (set of 2)
	191127	15-474-224	Adjustable foot rest. For ergonomic foot and posture support
	1910179	13-261-361	Shelf to hang from base stand
	1910083	13-261-362	Interior tubing retainer clips (6 per kit)
	ERG-CR-2330	ERG-CR-2330	Thermo Scientific Ergolign Saddle Stool, 23 inch to 30.5 inch adjustable seat height
<b>Sidewall Utilities</b>	1911308	15-474-208	Combustible gas tap for side wall
	1911309	15-474-200	Noncombustible gas tap for side wall
	1911310	15-474-210	Vacuum tap for side wall
	1911311	15-474-211	Water tap for side wall
	50059017	13-998-025	Hose connection to fit through side wall access port, fits 0.5 inch (10-13 mm) inner diameter tubing
	50076408	13-998-026	Hose connection to fit through side wall access port, fits 0.1 inch (3 mm) inner diameter tubing



50073700  
thimble duct exhaust transition



50078900  
utility basket



1911311  
water tap for side wall



Optional rear wall  
3 service tap connections

## factory-installed accessories

	Description
<b>Customized Models</b>	The Herasafe KS worksurfaces are fully customizable to accommodate any unique work task. Review ordering information with your Thermo Fisher Scientific sales representative. Lead times will be extended for custom orders.
<b>Worksurfaces</b>	KS 9 (3 ft) without worksurface — required for installation of non-standard worksurface options
	KS 12 (4 ft) without worksurface — required for installation of non-standard worksurface options
	KS 15 (5 ft) without worksurface — required for installation of non-standard worksurface options
	KS 18 (6 ft) without worksurface — required for installation of non-standard worksurface options
	Flat worksurface module 150 mm wide
	Flat worksurface module 300 mm wide
	Flat worksurface module 600 mm wide
	Flat worksurface module 900 mm wide
	Vibration damped worksurface 300 mm width. Constructed in stainless steel with granite plate. For weighing and microscopy equipment
	Vibration damped worksurface 600 mm width. Constructed in stainless steel with granite plate. For weighing and microscopy equipment
	Minisink tabletop insert 300 mm. Holds ~3 liters
	Indented worksurface module 600 mm wide
	Indented worksurface module 900 mm wide
	Indented worksurface module 1200 mm wide
Indented worksurface module 1500 mm wide	
Indented worksurface module 1800 mm wide	
<b>Universal Piping</b>	Rear wall 3 service tap connections, NPT 3/8, left side of cabinet
	Rear wall 3 service tap connections, NPT 3/8, right side of cabinet
	Rear wall 6 service tap connections, NPT 3/8, 3 on left and 3 on right side of cabinet. (May not be used with KS 9 (3 ft) cabinet.)
	Combustible gas tap for rear wall
	Noncombustible gas tap for rear wall
	Vacuum tap for rear wall
	Water tap for rear wall



50078875  
optional minisink tabletop insert



50076162 optional 600 mm vibration  
dampened granite block



Optional 316 stainless steel  
worksurface configurations


# Thermo Scientific 1300 Series A2

## specifications

Size / Width			3-foot	4-foot	5-foot	6-foot
<b>Dimensions</b>	Exterior dimensions WxHxD	inches	39.4 x 61.8 x 31.5	51.2 x 61.8 x 31.5	63 x 61.8 x 31.5	74.8 x 61.8 x 31.5
	Interior dimensions WxHxD	inches	35.4 x 30.7 x 24.8	47.2 x 30.7 x 24.8	59.1 x 30.7 x 24.8	70.9 x 30.7 x 24.8
	Shipping dimensions <sup>1</sup> WxHxD	inches	43.7 x 68.1 x 36.4	55.5 x 68.1 x 36.4	67.3 x 68.1 x 36.4	79.1 x 68.1 x 36.4
	Working height of front window	inches	Choice of 8 or 10			
	Maximum opening height of front window	inches	21			
<b>Weight</b>	Net weight	lbs	375	441	507	617
	Shipping weight <sup>1</sup>	lbs	419	496	573	694
<b>Electrical Data</b>	Voltage	V	120			
	Frequency	Hz	60			
<b>Filter Specification</b>			H14 HEPA EN 1822, 99.995% MPPS			
<b>Certification</b>			NSF/ANSI 49, UL, CE			
<b>Ergonomics</b>	Lighting power	fc	>100	>120	>120	>120
	SmartPort		Two 3" plugged cable ports, one on each side wall			
	Receptacles		2 GFI duplex			
	Service valves		Up to 6 total (3 on each side wall)			
<b>8-Inch Work Aperture</b>	Sound pressure level	dB (A)	<62	<62	<64	<64
	Energy consumption, operating set point <sup>2</sup>	W	150	180	280	360
	Energy consumption, Night Set-Back mode	W	70	70	120	120
	Heat output, operating mode (non-vented)	BTU/hr	512	614	955	1228
	Heat output, Night Set-Back mode (lights off, non-vented)	BTU/hr	239	239	410	410
	Exhaust/Inflow air volume	cfm	207	276	344	413
	Exhaust volume, thimble ducted	cfm	269	358	448	537
<b>10-Inch Work Aperture</b>	Sound pressure level	dB (A)	<63	<63	<65	<65
	Energy consumption, operating set point <sup>2</sup>	W	170	200	310	400
	Energy consumption, Night Set-Back mode	W	70	70	120	120
	Heat output, operating mode (non-vented)	BTU/hr	580	683	1058	1366
	Heat output, Night Set-Back mode (lights off, non-vented)	BTU/hr	239	239	410	410
	Exhaust/Inflow air volume	cfm	258	344	431	517
	Exhaust volume, thimble ducted	cfm	336	448	560	672

<sup>1</sup> Export packaging dimensions and weight upon request

<sup>2</sup> Clean filters, fans operating at nominal set point, interior lighting activated





## ordering info

Our 1300 Series A2 packages feature everything required to get started, including cabinet, manual adjustable height stand, factory-installed UV light and one set of armrests.

	Size	8-inch work aperture		10-inch work aperture	
		Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.
<b>SmartCoat Interior</b>	3-foot	1340	13-261-340	1321	13-261-332
	4-foot	1395	13-261-314	1385	13-261-304
	5-foot	1351	13-261-344	1369	13-261-336
	6-foot	1397	13-261-315	1387	13-261-306
<b>One-piece Stainless Steel Interior</b>	3-foot	1342	13-261-342	1323	13-261-334
	4-foot	1335	13-261-312	1375	13-261-222
	5-foot	1353	13-261-346	1371	13-261-338
	6-foot	1337	13-261-313	1377	13-261-223



Manual adjustable stand



UV light



Stainless steel armrests

To build your own 1300 Series A2, start with these cabinet-only part numbers:

	Size	8-inch work aperture		10-inch work aperture	
		Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.
<b>SmartCoat Interior</b>	3-foot	1339	13-261-339	1320	13-261-331
	4-foot	1325	13-261-310	1355	13-261-300
	5-foot	1350	13-261-343	1368	13-261-335
	6-foot	1327	13-261-311	1357	13-261-302
<b>One-piece Stainless Steel Interior</b>	3-foot	1341	13-261-341	1322	13-261-333
	4-foot	1305	13-261-308	1345	13-261-220
	5-foot	1352	13-261-345	1370	13-261-337
	6-foot	1307	13-261-309	1347	13-261-221

# Thermo Scientific 1300 Series A2

## customer-installed accessories

	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Description
<b>Stands</b>	1911400	13-261-347	Manual adjustable stand for 3 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	1911302	15-474-202	Manual adjustable stand for 4 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	1911401	13-261-348	Manual adjustable stand for 5 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	1911303	15-474-203	Manual adjustable stand for 6 foot cabinet, provides work surface height of 30 to 38 inches adjustable in 2 inch increments
	50119096	13-998-004	Electric adjustable stand for 3 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	50119126	15-453-300	Electric adjustable stand for 4 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	50119127	13-998-005	Electric adjustable stand for 5 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	50119128	15-453-301	Electric adjustable stand for 6 foot cabinet (120V), provides infinite height adjustability of 30 to 38 inches
	1911409	13-261-349	Rolling castor stand for 3 foot cabinet, provides fixed height of 30 inches
	1911322	15-474-325	Rolling castor stand for 4 foot cabinet, provides fixed height of 30 inches
	1911410	13-261-350	Rolling castor stand for 5 foot cabinet, provides fixed height of 30 inches
	1911346	15-474-326	Rolling castor stand for 6 foot cabinet, provides fixed height of 30 inches
<b>Exhaust Transitions</b>	1911316	15-474-216	Thimble duct exhaust transition for 3, 4 and 5 foot cabinet
	1911317	15-474-217	Thimble duct exhaust transition for 6 foot cabinet
	1910185	1910185	Alnor exhaust alarm for 1300 Series A2 biological safety cabinets, qualified install only
<b>User Convenience Equipment</b>	1910083	13-261-362	Interior tubing retainer clips (6 per kit)
	1910179	13-261-361	Hanging shelf for base stand
	1911313	15-474-213	Stainless steel armrest (set of two)
	1911127	15-474-224	Adjustable footrest
	1910112	13-261-363	Replacement SmartPort grommets (4 per kit)
	ERG-CR-2330	ERG-CR-2330	Thermo Scientific Ergolign Saddle Stool, 23 inch to 30.5 inch adjustable seat height
	1911413	13-261-351	IV bag holder kit with 6 hooks (for 3 foot cabinet only)
<b>Universal Piping Utilities</b>	1911312	15-474-212	IV bag holder kit with 12 hooks (for 4 foot, 5 foot, and 6 foot cabinet)
	249098	13-261-360	Stainless steel ball valve (replaces standard non-corrosive PVC drain valve)
<b>Sidewall Utilities</b>	1911308	15-474-208	Combustible gas valve
	1911309	15-474-209	Non-combustible gas valve
	1911310	15-474-210	Vacuum valve
	1911311	15-474-211	Water valve



1910179 hanging shelf for base stand



ERG-CR-2330 Ergolign Saddle Stool




191127 optional adjustable footrest

## factory-installed accessories

	Thermo Scientific Cat. No.	Fisher Scientific Cat. No.	Description
<b>Light Option</b>	1911319	15-474-219	UV light option for 3 foot or 4 foot cabinet
	1911368	15-474-327	UV light option for 5 foot or 6 foot cabinet
<b>Universal Piping</b>	1910286	13-261-352	For 3 foot <b>coated</b> interior models, right side
	1910287	13-261-353	For 3 foot <b>coated</b> interior models, left side
	1911321	15-474-221	For 4 foot <b>coated</b> interior models, right side
	1911369	15-474-322	For 4 foot <b>coated</b> interior models, left side
	1910288	13-261-354	For 5 foot <b>coated</b> interior models, right side
	1910289	13-261-355	For 5 foot <b>coated</b> interior models, left side
	1911367	15-474-323	For 6 foot <b>coated</b> interior models, right side
	1911370	15-474-324	For 6 foot <b>coated</b> interior models, left side
	1910290	13-261-356	For 3 foot <b>stainless</b> interior models, right side
	1910291	13-261-357	For 3 foot <b>stainless</b> interior models, left side
	1911381	13-261-224	For 4 foot <b>stainless</b> interior models, right side
	1911382	13-261-225	For 4 foot <b>stainless</b> interior models, left side
	1910293	13-261-358	For 5 foot <b>stainless</b> interior models, right side
	1910294	13-261-359	For 5 foot <b>stainless</b> interior models, left side
	1911383	13-261-227	For 6 foot <b>stainless</b> interior models, right side
	1911384	13-261-228	For 6 foot <b>stainless</b> interior models, left side
	<b>Universal Piping Utilities</b>	1911304	15-474-204
1911305		15-474-205	Non-combustible gas valve for universal piping option
1911306		15-474-206	Vacuum valve for universal piping option
1911307		15-474-207	Water valve for universal piping option

Provides two pre-plumbed ports with connections to top, side or bottom of the cabinet



Universal piping

# Your choice for containment, comfort and convenience



Quick Comparison Chart		Traditional BSC	Herasafe KS	1300 Series A2
<b>Containment</b>	SmartFlow dual motor design		•	•
	Digital Airflow Verification (DAVe) monitoring		•	•
	Night Set-Back mode (maintains sterility)		•	•
	H14 HEPA filters (99.995% efficiency at MPPS)		•	•
	Negative pressure side wall plenums		•	•
	All metal plenums		•	•
	Performance Factor/SmartFlow Indicator		•	•
	One-piece indented work tray		•	
	Cross-beam UV disinfection system		•	
SampleGuard aerosol-tight sealing window frame		•		
<b>Comfort</b>	Low noise level		•	•
	10° sloped front		•	•
	Bright workspace		•	•
	Side glass walls		•	
<b>Convenience</b>	Programmable UV light		•	•
	Eye-level control panel/display		•	
	SmartClean Plus hinging window		•	
	SmartClean window cleaning			•
	SmartPort tubing/cable organizer			•
	Motorized front sash		•	
	316 stainless steel interior		•	
<b>Efficiency</b>	Energy efficient DC motors		•	•
	Night Set-Back mode (reduces energy)		•	•
	Extended motor life (>100,000 hours)		•	•
	8 inch window option			•



© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



**In the United States:**

For customer service, call 1-800-640-0640

To fax an order, use 1-800-290-0290

To order online: [www.fisherhealthcare.com](http://www.fisherhealthcare.com)

BN0812161

BRBSC1300KSHMD 0816

