



Thermo Scientific SPECTRONIC 200
Visible Spectrophotometer



The perfect
teaching instrument



Designed for the Teaching Laboratory

Thermo Scientific™ SPECTRONIC™ spectrophotometers have served as core analytical instruments in school, college and university teaching laboratories since 1953. The SPECTRONIC 200 visible spectrophotometer continues this tradition as the new standard for the 21st century. Robust, innovative features ensure an excellent teaching instrument that will last. Delivering quick, convenient measurements in the hands of novice users, the SPECTRONIC 200 spectrophotometer is the perfect fit for the teaching lab.

Smart Technology, Simplicity and Speed

Smart Technology

The CCD detector measures each portion of the spectrum with the ideal integration time to account for lamp intensity and detector sensitivity. This provides the best possible photometric performance for accurate, reliable results no matter where in the spectrum a measurement is made.

Simplicity

- No need to do a zero %T measurement – the instrument does this at start-up
- No need to repeat a 100%T measurement every time you change the wavelength in Live Display mode. When you *autozero*, 100%T is recorded at all wavelengths
- Wavelength knob labeled with a λ to set or change the wavelength
- One touch printing with a USB printer

Speed

- Live measurements update every 2 seconds
- Scan from 400 nm to 900 nm in about 10 seconds

Total Control without a Computer

The stand-alone SPECTRONIC 200 spectrophotometer provides the experiment modes you use with its own controls and screen that offer:


- Live display with color spectrum indicator
- Full spectrum scanning for peak identification
- Emulation of legacy instrument interfaces for easy integration into existing laboratories and protocols
- Quantitative analysis with up to four standards or a user entered factor
- Wavelength selection with the λ knob or the arrow keys
- Four point fixed wavelength measurements

Adjustable Tilt Screen

Set the screen at the perfect angle for viewing by users of all heights. Fold it down below the top level of the instrument for compact and safe storage.


Coarse/Fine Wavelength Setting

Turn the λ knob to change the wavelength in steps of 10 nm. Press it down and turn to change the step to 1 nm. The range is from 340 nm to 1000 nm, and there's no hard stop for the knob at the ends of the range.




Classroom Friendly Sample Compartment

Whether you measure in 10 mm square cuvettes or in test tube cuvettes up to 25 mm diameter, the standard sample compartment stage adapts to be the perfect fit for your sample.




Students spill chemicals. You can't prevent it, but the SPECTRONIC 200 spectrophotometer is ready to handle it. The sample compartment features:

- Capacity for up to 250 mL of spilled liquid
- Simple lift-out for easy cleaning in the sink
- Plastic construction to resist acid solutions
- Lid that covers the control pad and screen for extra protection when samples are being moved in and out of the compartment



Sealed Keys

Tactile rubber keypad forms a seal with the top cover to prevent spills from entering the instrument.





Sloping Surfaces

There are no horizontal surfaces on the SPECTRONIC 200 spectrophotometer that invite students to balance their samples on it. All surfaces on the top cover slope to shed spilled liquids.

Convenient Cuvette Rack Compartments

Compartments to the right and left of the sample compartment are equipped with removable cuvette racks designed to hold six square or test-tube cuvettes. Have your students prepare their samples in the wet chemistry area and put them in a rack. Place the full rack in one compartment and the empty rack in the other. It's easy to keep clear track of which sample is which as they measure and move samples from one side to the other.




Modern On-Board Software



Live Display for the 21st Century

Measurement wavelength and result are shown in large clear text. A spectrum along the bottom of the screen indicates what color this wavelength represents. The display updates regularly, and students can freeze the display if desired.



Scan

Scan a full spectrum or a defined range. Use the cursor to see precise absorbance values.

Multi-wavelength


Perform advanced experiments where you measure multiple wavelengths simultaneously.

λ 1	430 nm	0.17	ABS
λ 2	610 nm	0.13	ABS
λ 3	339 nm	-----	ABS
λ 4	339 nm	-----	ABS

Enter to HOLD

Quantitative Analysis with Standards or a Factor


Make a Beer's Law plot with up to four standards and print it at the touch of a button. The SPECTRONIC 200 spectrophotometer shows your students their result and where their sample lies on the calibration curve.



UV-VISIBLE SPECTROPHOTOMETRY





Expanded Capabilities



Full Power of the Instrument at the Click of a Mouse

Thermo Scientific VISIONlite™ software offers a modern look with intuitive controls that students will find easy to use. They'll be off and running VISIONlite software with minimal instruction.


The experiment parameters are clearly displayed, offering only the settings you need. There are no confusing options to distract students or spoil data.



Data display is clear and concise. Graphical display includes optional peak picking and curve fitting. Data tables are displayed below the graph in a scrolling frame.



Quantitative analysis gives you flexibility to fit multiple curve types or enter a factor. Customize your data display precision and result unit. Result reports include the calibration graph, curve fit, standard solution data and experimental results.



Use Your Existing Protocols

Emulation Modes Offer Seamless Integration


Whether you are equipping a whole laboratory or adding to your existing instrument set, the SPECTRONIC 200 spectrophotometer smooths the transition from old to new technology. The SPECTRONIC 200 on-board software includes full emulation of SPECTRONIC 20 and Thermo Scientific GENESYS 20 control sets. This allows you to continue using your current laboratory protocols and instruction sets. Plus, if budgets are tight, you don't need to replace all your instruments at once.

SPECTRONIC 20 Emulation

All data modes from this Thermo Scientific stalwart are faithfully reproduced in this emulation. Right and left arrows change mode, the knob sets wavelength, and the up/down arrows set concentration and factor for quant methods. It's just like using a SPEC 20 instrument.

GENESYS 20 Emulation

Use the navigation pad to highlight and press on-screen buttons. Your existing operating instructions will work perfectly for students running this emulation mode.



The Perfect Fit – On or Off the Bench

The SPECTRONIC 200 spectrophotometer takes design for the teaching laboratory to a new level. We've put the experience and knowledge gained from supporting over half a million SPECTRONIC instruments into designing an instrument that's worthy of the SPECTRONIC name and ready to serve in classrooms around the world. The SPECTRONIC 200 spectrophotometer makes teaching easy in so many ways:

Easy to Set up

- Just plug it in and power it up
- No computer or data logger to connect and maintain

Easy to Use

- Supports cuvettes or test-tubes up to 25 mm in basic configuration
- Room-light resistant – run tall test tubes with the lid open
- Full color, variable angle LCD screen
- New scan and multi-point quant software
- Legacy instrument emulation modes
- Intuitive navigation
- Convenient coarse/fine wavelength control
- Includes cuvette racks and storage compartments

Easy to Maintain

- Sample compartment lid protects controls from spills while open
- Sample compartment contains up to 250 mL of spilled liquid and lifts out for easy cleaning

- All top surfaces slope to shed spills
- Polymer construction for chemical resistance
- No motorized parts to wear out
- Fixed grating for consistent wavelength accuracy
- Trap-door mounted lamp for easy replacement

Easy to Store

- Rubber feet on the back surface let you stack instruments on a shelf like books
- Storage compartments hold your empty cuvette racks
- Removable power cable can be coiled and stored in the sample compartment





Specifications

Optical Design		Single Beam
Spectral Bandwidth		≤4 nm
Light Source		Tungsten-halogen
Detector		2048 element CCD
Wavelength	Range	340 nm to 1000 nm
	Accuracy	±2 nm
	Repeatability	±1 nm
	Data interval	1 nm
Photometric	Range	-0.3 A to 2.5 A
	Readout	ABS, %T, Concentration
	Accuracy	±0.01 A at 0.3 A ±0.05 A at 1.0 A (SPECTRONIC standard filters measured at 590 nm)
	Repeatability	±0.3 %T at 50 %T
Stray Light		<0.2 %T (with SPECTRONIC standard SRM 400 filter)
Display		Variable angle 320 × 240 pixel color graphical LCD 7 × 5 cm, 8.6 cm diagonal (2.75" × 2", 3.4" diagonal)
Keypad		Sealed tactile rubber
Standard Features	Sample compartment	Lifts out for cleaning
	Cuvette racks and compartments	2 included with dedicated compartments
	Included cuvettes	1 cm plastic (quantity 10)
Standard Interfaces		USB-B for connection to a remote computer USB-A for connection to a printer or USB memory device
Languages		On-board software in English, Spanish, French, German and Italian
Power Requirements		100–240 V, 50–60 Hz (selected automatically)
Dimensions		39 cm W × 30 cm D × 16 cm H (15.3" W × 11.8" D × 3.6" H)
Weight		4.4 kg (9.7 lbs)
Warranty		1 year



Ordering Information

Description	Part Number
SPECTRONIC 200, US Power Cord	14-385-491
Optional Software and Accessories	Part Number
VISIONlite 5 Software, Single User License	14-385-500
VISIONlite 5 Software, Six User License	14-385-501
VISIONlite 5 Software, Twelve User License	14-385-502
Long Path Rectangular Cell Holder to 100 mm	14-385-493
Long Path Cylindrical Cell Holder to 100 mm	14-385-494
Cuvette Rack – package of 6	14-385-495
Epson® TM-T88 Continuous Feed Printer	840-234000
Replacement Lamp	14-385-499



Part of Thermo Fisher Scientific

For customer service, call 1-800-766-7000.
To fax an order, use 1-800-926-1166.
To order online: www.fishersci.com

BR52003_E 01/13M BN1128126



Part of Thermo Fisher Scientific