

biotechne®

R&D SYSTEMS

# DuoSet™ ELISA Development Systems

## Limitless Possibilities

### Customizable to Your Needs

Our development systems are adaptable to multiple platforms, including Gyrolab™, MSD®, Luminex®, and 384-well plates. We also offer tailored protocols and ancillary kits.

### One-stop Shopping

We offer the largest selection of analytes (>1000) for the largest variety of species.

### Trustworthy Performance

We are the most referenced manufacturer of immunoassays because of our unparalleled performance and reliability.

### Consistent Supply

We can guarantee our supply since we manufacture all reagents in-house.


### Economical

5-plate packs are available for many analytes.



## Trustworthy Performance

### Nearly identical performance—measured 10 years apart



Comparison of Human  $\beta$ -NGF Standard Curve O.D.s from 2003 and 2013. Using the Human  $\beta$ -NGF DuoSet ELISA Development System (Catalog # DY256), standard curve values generated in 2003 and 2013 were compared for lot-to-lot consistency. Standard curve O.D.s remained consistent over ten years.

### DuoSet ELISA Kit Components

- ELISA capture reagent
- Biotinylated detection antibody
- Streptavidin-HRP
- Mass-calibrated standard
- Detailed protocol

### Ancillary Kit Components

- 96-well plates
- Plate sealers
- Substrate reagents
- Stop solution
- ELISA coating buffer
- Reagent diluent
- Wash buffer

Distributed by Fisher Scientific. Contact us today:

In the United States

Order online: [fishersci.com](http://fishersci.com)

Call customer service: 1-800-766-7000

© 2022 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at [fishersci.com/trademarks](http://fishersci.com/trademarks).

BN224636552-US

 **fisher scientific**  
part of Thermo Fisher Scientific